

Protocol: BCA Total Protein Assay

Notes:¹

- Protocol was modified from: Thermo Scientific Protocol and Nanodrop 1000 BCA protocol
- See additional Thermo Scientific BCA kit information appended to protocol for buffer compatibility, nanodrop 1000 capability, and kit sensitivity.

Materials:

1. 1x PBS, prepared by N/A on N/A, location: stored at room temperature in 1-L glass reagent bottles rm. 231
2. Lysis Buffer, 300 mM NaCl, 50 mM Tris, pH 7.4, 0.5% Triton X-100, Roche anti-protease cocktail, location: stored at room temperature rm. 231
 - a. Sodium Chloride, Fisher Scientific, Lot# 094963, CAS# 7647-14-5, received 2/10/10, location: reagent shelf rm. 231
 - b. Trizma Hydrochloride, Sigma, Lot# 78H5403, Received: 1/28/99, opened: 10/26/99, location: reagent shelf rm. 231
 - c. Triton X-100, Sigma, Lot# 52H0286, received 7/93, location: reagent shelf rm. 231
 - d. Roche anti-protease cocktail (currently not available in lab)
3. Thermo Scientific BCA Protein Assay Kit, Prod# 23227, Lot# JJ126728, received N/A, opened N/A, location: reagent shelf rm. 231
4. PCR thermocycler, location: rm. 211
5. PCR tubes, Simport, Lot# 12211865
6. Nanodrop 1000
7. Ice

Methods:

BCA Protein Assay

1. Prepare BSA standards using PBS as a diluent
 - a. Standard curve range for 1:1 ratio assay is approximately between 10 µg/mL and 200 µg/mL.
 - b. Standard curve range for 1:20 ratio assay is approximately between 200 µg/mL and 8000 µg/mL

¹ Authors/Revision: Created by SP on 2012-07-17; Modified by SP on 2012-08-12

- c. A “0” standard/reference must be made using PBS in place of sample
 - d. Selection of standard curve range depends on the source of the sample
- 2. Prepare working reagent
 - a. 50 parts reagent A + 1 part reagent B
 - b. Amount of working reagent per tube is 10 μ L/tube for 1:1 ratio assay, and 80 μ L/tube for 1:20 ratio assay
- 3. Aliquot proper volume of standards and samples into PCR tubes
 - a. Volume of standard/sample per tube is 10 μ L/tube for 1:1 ratio assay and 4 μ L/tube for 1:20 assay
- 4. Add proper volume of working reagent to each tube
 - a. Volume of working reagent per tube is 10 μ L/tube for 1:1 ratio assay and 80 μ L/tube for 1:20 assay
- 5. Incubate tubes in PCR thermocycler at 60 degrees C for 30 minutes.
 - a. Program is named “BCA”.
 - b. If evaporation is a problem, incubate tubes at 37 degrees C for 30 minutes.
See additional kit information from Thermo Scientific for full details.
- 6. Quantify protein concentration using Nanodrop 1000 BCA program
 - a. Use 2 μ L of sample/standard on Nanodrop sensor
 - b. Use water as the blank
 - c. Reference measurement is the “0” standard
 - d. Remember to save standard curve and data in Gregson folder.