

Flow Microwell Staining—Dendritic Cell Panels

Supplies:

1. 200 mcL round-bottom plate
2. PDBS + 4% FCS (washing solution)
3. PFA 1% (poison)
4. antibodies
5. microcentrifuge tubes to mix antibodies
6. flow tubes with lids

Table 3: Dendritic Panel #3. March 20, 2008 onward.

FITC	FL1-H	BDCA-2 / CD303
PE	FL2-H	BDCA-1 / CD1c
PerCP	FL3.H	CD45

Table 4: Dendritic Panel #4. March 20, 2008 onward.

FITC	FL1-H	HLA-DR
PE	FL2-H	CD14
PerCP	FL3.H	CD45

Panels

Table 1: Dendritic Cell Panel #1: Used March 14/18/19.
Was panel #2 on 3/14.

FITC	FL1-H	BDCA-2 / CD303
PE	FL2-H	CD14
PECy5	FL3.H	HLA-DR
APC	FL5 (FL2-W)	BDCA-4 / CD304

Table 5: Dendritic Panel #5. March 20, 2008 onward.

FITC	FL1-H	HLA-DR
PE	FL2-H	CD14
PerCP	FL3.H	CD45
APC	FL5 (FL2-W)	BDCA-3 / CD141

Table 2: Dendritic Panel #2. March 20, 2008 onward.

FITC	FL1-H	BDCA-2 / CD303
PE	FL2-H	CD14
PerCP	FL3.H	CD45

Table 6: Dendritic Panel #6. March 20, 2008 onward.

FITC	FL1-H	HLA-DR
PE	FL2-H	BDCA-1 / CD1c
PECy5	FL3.H	CD45
APC	FL5 (FL2-W)	BDCA-3 / CD141

Surface Staining:

(*< 60 minutes*)

1. Get sample into each well (each well will become a flow tube)
 - 1.1. place 200 mcL of sample into each well (dilute with washing solution)
 - 1.2. spin plate 1200rpm x 2 minutes
 - 1.3. flip drip dry
 - 1.4. wipe off top
 - 1.5. agitate via quick vortex
2. Stain the cells
 - 2.1. add 16 mcL of washing solution
 - 2.2. add 2.5 mcL of each antibody to wells
 - 2.3. mix with 10 mcL pipet (up and down 5—10 times)

- 2.4. incubate at 4 °C for 30 minutes
3. Wash antibody off of cells
 - 3.1. add 200 mcL of washing solution
 - 3.2. spin plate 1200rpm x 2 minutes
 - 3.3. flip drip dry
 - 3.4. wipe off top
 - 3.5. agitate via quick vortex
4. Fix cells
 - 4.1. add 200 mcL 1% PFA to each well
 - 4.2. fill each flow tube with 300 mcL PFA
 - 4.3. add contents of wells to tube
 - 4.4. vortex tubes
 - 4.5. cover with cap, keep at 4 °C until ready to use